
Benex®
Pat.No. CH 696 458

www.benex-dent.com

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Anaesthesia. Cutting periodontal fibres (Sharpey fibres) in the sulcus by using the **ZEPP** Periotome **Xiod**.

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Strong, large roots must be loosened/luxated by axial movements within 30 sec., using a slim elevator/twister (**Xiod** from **ZEPP**). Without using transversal movements. The roots of multiradicular teeth have to be separated and extracted individually.

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Drilling with the diamond coated twist drill (choose the drill according to the used screw) should be in the axis and center of the root fragment. It should be approx. 7 mm in the hard tissue, deeper drilling will not be necessary. Drilling is performed with water-cooling. In order to remove drilling chips more easily, an inward and outward movement is recommended for deep drilling.

Recommended rpm: 500 - 700 rpm. | 12.300.30 max. 3000 rpm. | 12.300.35 max. 2200 rpm.

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According to circumstances, the short or long (thin or thick) extraction screw with screwing support is inserted.

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The Extractor is positioned on the adjacent crowns: The opening of the round, revolvable segment plate is adjusted in vestibular direction ensuring a good view of the extraction screw. After the traction rope has been hooked into the extraction screw, it is guided over the reverse roller and fixed to the hook of the extraction slide. Under slight traction – so that the rope does not hang out – the instrument is placed on the adjacent teeth by turning the hand screw. During positioning it is important to see that both the screw and the rope do have the same axial direction.

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Once the extractor is positioned properly, the extraction is carried out by turning the hand screw. In the case of strong, long roots the periodontal fibers have to be prestressed with submaximal tensile force during 30 - 40 seconds. This means that the instrument is positioned in such a way that it has a support under tension. Then you have to wait in this position. Within 30 - 40 seconds the fibers will yield, and thus the instrument can be prestressed again. This procedure has to be repeated, if necessary.

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The Free-End Situation: In order to guarantee an optimized, vertical power-transmission, a right-angle level has to be available in order to apply the Benex® Extractor. With the new developed quadrant support the supporting level can be corrected and/or a bad level can be amended. This is often necessary with Benex® extractions in the mandible and with palatinal roots in the maxilla.

1. Fill the quadrant support with a fast-setting plasticine and model a gap in the buccal hole region.
2. Place the quadrant support with opening in the root extraction area, at right angles with the root axis.
3. Remove the quadrant support when the plasticine begins to set. The complete setting happens extra-orally.
4. Root luxation. Afterwards insertion of the Benex® Screw and repositioning of the quadrant support and the Benex® Extractor.

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

In order to achieve an optimal cleaning result, the Benex-Extractor has to be completely disassembled before subjecting it to the cleaning process. The disassembly can be carried out as follows:

First of all, remove the 'knurled screw' from the extractor by turning the knurled screw anti-clockwise off the thread.

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
In a second step the 'swivel segment' is removed from the base plate by turning it anti-clockwise or clockwise.

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
Step 3 is to remove the 'cross pins' from the base plate in a slight pulling or pushing movement.

IMPORTANT: Please fix the guide cover of the threaded pin securely during this procedure in order to prevent the integrated pressure spring from unintentional extension.

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
 In step 4, please remove the guide cover [1] with care. In doing so, fix the ‘tensioning element’ [2]. Remove the pressure spring [3] first and the ‘tensioning element’ [2] subsequently. Also remove the ‘O-ring’ gaskets [4] from the holders. During the assembly these have to be inserted prior to the cross pins in order to guarantee a precise holding.

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 Step 5 is to remove the ‘POM plate’ from the base plate by pulling slightly.

Afterwards, check the correct fit of the ‘O-ring’ gaskets.

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 The last step is to remove the ‘POM segment plate’ from the ‘swivel segment’ by slight pull.

Thus, the components can be cleaned according to requirements and be assembled afterwards.

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