Handelsname: ZEPF Kronenspreizzange REF 19.277.01Z

Trade Name: ZEPF Crown-Spreading Pliers REF 19.277.01Z

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Dear Customer!

We thank you very much for having chosen this high quality instrument from Helmut Zepf.

To ensure that you are satisfied with the instrument for a long period, please find a few instructions for use below.

- 1) First separate the crown completely using the crown-spreading pliers. Please ensure that the crown has been split deep enough. Sharply defined edges with the separator (as cylindrical as possible) give you secure adaption in the slot and greatly facilitate subsequent spreading.
- 2) Now the tip of the crown-spreading pliers is inserted into the slot from occlusal, buccal or lingual, depending on the access option.
- 3) The good anatomical curvature of our crown-spreading pliers allows their use in all quadrants and thus guarantees safe and correct use.
- 4) Careful, directed and slow spreading of the crown now allows widening of the separating gap. This should be at least 1 mm. Compared to normal crown removers, do not apply leverage! This way no pressure is applied to the root, the adjacent tooth or the jaw. The point of rotation is at the tip of the instrument and protects the root underneath.
- 5) Then the crown can be removed without difficulty.

An Allen wrench is provided to position the tips in parallel. If the screw becomes loose during cleaning, the working end can be repositioned after cleaning and before sterilization.

Please note that the wrench is not intended for autoclaving.

Download 
ZEPF|care
 Optimales
 Reinigen

Download 
ZEPF|care
 Optimum
 Cleaning