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Benex[®]
PENDING PATENT

🇩🇪 FAQ | 🇬🇧 FAQ

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
The Alveolar Wall after Benex[®] Extraction

Documentation acc. to Benno Syfrig: In my opinion, the protection / conservation of the buccal alveolar wall by longitudinal extraction is very advantageous and important. Apparently, the micro structures of hard / soft tissues are less injured and the osseous ridge is conserved, respectively osteolytic parts are reossificated.

In the past two weeks, I took photos of this frequently observed reossification / recontouring of the alveolar wall. In the first six to twelve weeks after extraction (depending on the size of the alveolar wall defect), newly built bone is sticking to the perioste when opening. Later, you will find a reossificated ridge and a smooth perioste.

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
Case 1 / Picture 1: Situation eight weeks after extraction. During extraction, a five millimetres deep buccal alveolar wall defect was photographed. The reossification of the buccal alveolar wall and bone fragments sticking to the perioste are well visible, see picture 1.

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Case 2 / Picture 2 and Case 3 / Picture 3 and 4: Three months and more after extraction, the smooth perioste can be relieved from the reossificated defect. In both cases I found - as in case 1 of the Benex[®] Extraction - a five millimetres deep, buccal alveolar wall osteolysis. Due to a patient delay the implantation was effected more than three months later. In both cases a smooth perioste was found. In case 2, picture 2, the extension of the alveolar wall reossification can be well seen.

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
Case 4 / Picture 5 and 6: In the following patient, an intact buccal alveolar wall was found during extraction. In this case, no periosteal reossification was found or was necessary. Two months after extraction an intact alveolar wall was found and the perioste was smooth.

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Case 5 / Picture 7: The following patient is a negative example: The referring dentist extracted gently! Three months after extraction a big alveolar wall defect was found. There are no signs of reossification on the perioste.

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The Free-End Situation

In order to guarantee an optimized, vertical power-transmission, a right-angle level has to be available in order to apply the Benex® Extractor. With the new developed quadrant support the supporting level can be corrected and/or a bad level can be amended. This is often necessary with Benex® extractions in the mandible and with palatal roots in the maxilla. The first 3-4 times the Benex shall be inserted in the front with unproblematic support for teeth. The quadrant support with plasticine serves as supporting area/counter bearing for extraction in the free-end gap. Due to the elastic counter bearing (plasticine/mucous membrane), the extraction force is slightly reduced. For this reason, the root to be removed has to be luxated more than with a usual Benex Extraction with support for teeth.

The supporting level for the Benex has to be at right angles with the axis of the teeth-/screw in order to apply vertical force. With the quadrant support the supporting level for the Benex can be corrected and/or a bad level can be improved. This is often necessary with Benex extractions in the mandible and with palatal roots in the maxilla.

Using hints:



1. Fill the spoon with quick binding plasticine, press the plasticine away with a finger within the areas of the buccal spoon opening.
2. Take the spoon into position: spoon opening in the root extraction area, at right angles with the root axis.
3. Remove the spoon from the mouth when the plasticine starts to bind. Finished bind, extra oral.
4. Root luxation / Screw assembly.
5. Repositioning of the Benex spoon, Positioning of the Extractor.

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The highly carious root

Before the Benex drilling, the carious, soft dentin should be removed with a bud bur. The cross section of the root fragments should be displayed, that the drilling stops in the centre of the fragment.

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
The twisted root

The root must be luxated well intraalveolarly.

The single rooted tooth: Quite frequently, the curved root breaks off during the Benex Extraction in the area of the strongest curve or on the height of the drill tip. In the second Benex drilling procedure, or if the curve of the root is very strong, after osteotomy the apical fragment can be turned/extracted as per the radius of the curve, using the root extractor. The osteotomy in the apical root part is not necessary for the reossification of the alveole, respectively for the preservation of the alveolar ridge.

For multirooted teeth it is recommended to extract the straight roots with the Benex first. This makes way for the conventional luxation of the curved root, in accordance with the radius of the curve. The palatinal roots of the maxillar molars are usually straight and can be extracted with the Benex, and the buccal roots can be rotated/extracted with the lever. For the Benex extraction of mandibular molars, the distal root (usually the tap-root) is more suitable. The wide, often small, sand glass shaped mesial root can then be bent/luxated conventionally in distal direction. In the majority of cases, this root is not suitable for the Benex screw fixation.

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The Germectomy

I understand that the dentists have doubts with the Benex® extraction on the germectomy. It did not exist a system until today which allowed a bore in the dental enamel and allowed there a self cutting screw fixation. This is actually new and the reason why the bur and the screw design have been patented as well.

The bur is so sharp that in the flexible tooth germ (swimming in the follicle) the bur hole can be actually made without problems and the self cutting screw can be turned in. As a matter of course the hole through the bone must be big enough for the place where the crown comes on. Either it has to be made an osteotomie or the crown must be split, whereas the bur hole is made advantageous for the screw. Take the long screw without pull rope-instrument and pull by hand the screw, the tooth germ, out of the follicle. This handling is also essential for the mesiodens.

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Problems with Pilot Drilling

One problem that I often meet, concerns the pilot-/inlay bur. Especially the cutting output of the diamond tips is insufficient. I need relatively much pressure to enforce the pilot drilling rather efficiently. The preparation of the inlay canal is also affected. I already changed to Endodontic Instruments, because the processing is more efficient and faster. Is there anything I am doing wrong or am I not the only one to face this problem? I use this system many times, but not as often as requested, please note problem above.

Suggested solution from Benno Syfrig, Dr. med., med. dent.

This problem is actually known. During drilling in a root treated canal the drill surface often daubs and the drill performance drops off. Therefore I recommend to remove the root filling material with an Endo-Instrument and then to use the Benex drill. With apical root rests where the canal is hard to find, this is not necessary. The canal is small there and the drilling is within the dentin. Furthermore the problem can be solved by means of the thicker screw. This screw does not have to be turned in that deep.

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When screw pivoting the root fragment splits

Possible causes and counteractive measures:

1. The drilling is not in the fragment centre. In order to hit the fragment centre, it is recommended to show the cross section of the fragment by means of a bud bur, resp. to remove soft tissues prior to the Benex drilling

2. The screw thread is blunt. For the extraction of delicate root fragments, the screw thread has to be sharp, as a blunt screw thread has a higher explosive force (see topic „Screw pulls out“).

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Screw pulls out!

Possible causes and counteractive measures:

1. Drilling/Screw lays in the widely prepared root canal. Use the thicker Benex screw or put the thin, long screw apically deeper into hard dentin. For this purpose, the drill extension might be required

2. Drilling/Screw has a perforated fragment wall or lays in soft, carious dentin. When drilling, you feel the moment of the perforation, the resistance gives up. Make a new drill apically deeper and if possible in the centre of the fragment. To hit the centre it is recommended to remove the fragment substance (carious tissue or root canal filling material) with a bud bur prior to the Benex drilling.

3. The screw thread is blunt. One Screw can be used more than twenty times. Check the sharpness of the thread from time to time: clench the screw thread between thumb and index finger and pull the screw thread with the thumb and index finger of the other hand. A screw with good working thread gets caught. A screw with a blunt thread is missing the retention, it has to be replaced.

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
Root does not separate!

Possible causes and counteractive measures:

1. The pull rope is not adjusted according to the root axis. Correct the extracting position (pull rope according to the root axis), at best by means of the quadrant support.

2. The root is insufficiently luxated intraalveolarly. You can leave the extractor under high efficacy and press the root in mesial / distal direction with a lever.

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