















**KlingEX<sup>REF</sup> 46.005.70**


-  **Sicherheitsklingenabnehmer**
-  **Safety Blade Remover**
-  **Changeur de lames de bistouri**
-  **Remover hojas de bisturís seguramente**
-  **Leva lama di sicurezza**



-  **positionieren**
-  **positioning**
-  **positionner**
-  **posicionar**
-  **posizionare**



-  **Klinge anheben durch Pressen**
-  **lift blade by pressing**
-  **élever la lame en pressant**
-  **elevant la hoja por presionar**
-  **lama levare tramite pressione**



-  **abziehen**
-  **pull off**
-  **retirer**
-  **descontar**
-  **estrarre**