

Pi
SERIES


Prosthetic Kit



Pi
SERIES

Distributed by:

Produced by:




HELMUT ZEPF
MEDIZINTECHNIK GMBH

78606 Seitingen-Oberflacht, Germany | Tel.: +49 (0) 74 64 / 98 88 -0 | Fax: +49 (0) 74 64 / 98 88 -88
E-Mail: info@zopf-dental.com | Internet: www.zopf-dental.com


HELMUT ZEPF
MEDIZINTECHNIK GMBH
www.zopf-dental.com

Prosthetic Kit

47.830.00 Complete Kit



The **ZEPP** Prosthetic Kit will simplify your work essentially! This is applicable for more than 90% of the existing abutments on the market. Standardly the kit will be delivered with a ratchet without torque function. With a click-function the handle is easily demountable. All inserts can be picked up with a contra-angle handpiece with Hex adapter (RA-HEX). A ratchet with torque function is available optionally. REF 47.803.02

The instrument rack is compatible to our Tray System (Catalogue Page 10-1 to 10-03) and can be used universally in all Washbasket sizes. Thus a platform is created and there is no limit to combine and enlarge it with additional instruments and tools.

A safe transport / accommodation according to the accident prevention regulation 103 and the cleaning solution as per RKI-Guidelines is consequently given.

The RA-Hex connection facilitates the adaption of the inserts in a medical handpiece.

Based on the RA-Hex connection is the Universal-Bone-Condensing-Kit, ref. 47.937.00 additionally available. It originates an added value under consideration of the used ratchets, insertion tools etc.



47.525.50
Driver Guide

47.525.55 Ratchet with demountable handle

Optional 47.803.02 Torque Wrench

Prosthetic Kit



Article description

47.830.00 Prosthetic Kit consisting of:

85.182.50 1/3 Washtray System

47.830.01 Rack for Prosthetic Kit

Order quantity

[]

[]

[]

Prosthetic Kit

Advantages



In the RA-Hex-Adapter the inserts will be picked up.












RA-Hex Adapter inserted in ratchet. Shown with demounted handle for use as finger-ratchet. Information: The RA-Hex-Adapter have a rotating finger rest available (see red arrow).















Among the recess for two RA-Hex-Adapters a socket for a comfortable removal of the screwdriver-inserts from the RA-Hex-Adapter is located.

▼ SHORT Screwdrivers with dental lock

	47.832.01 TORX T6, Straumann Aesthura	[]
	47.832.02 Universal, flat 1.6 mm	[]
	47.832.03 Universal, flat 2.0 mm	[]
	47.832.04 HEX 0.03", Allen Key, Camlog	[]
	47.832.05 HEX 0.05", Allen Key, Camlog, Sulzer (Zimmer), Semados, 3I	[]
	47.832.06 HEX A/F 0.9 mm, Allen Key, IMPLA, TIOLOX, BREDDENT, XIVE	[]
	47.832.07 HEX A/F 1.0 mm, Allen Key, Ankylos	[]
	47.832.08 HEX A/F 1.2 mm, Allen Key, IMPLA Nobel Biocare, Frialit, XIVE, IMZ, Biomet 3I	[]
	47.832.09 HEX A/F 1.8 mm, Allen Key, Ankylos	[]

▼ LONG

	47.833.01 TORX T6, Straumann Aesthura	[]
	47.833.02 Universal, flat 1.6 mm	[]
	47.833.03 Universal, flat 2.0 mm	[]
	47.833.04 HEX 0.03", Allen Key, Camlog	[]
	47.833.05 HEX 0.05", Allen Key, Camlog, Sulzer (Zimmer), Semados, 3I	[]
	47.833.06 HEX A/F 0.9 mm, Allen Key, IMPLA, TIOLOX, BREDDENT, XIVE	[]
	47.833.07 HEX A/F 1.0 mm, Allen Key, Ankylos	[]
	47.833.08 HEX A/F 1.2 mm, Allen Key, IMPLA Nobel Biocare, Frialit, XIVE, IMZ, Biomet 3I	[]
	47.833.09 HEX A/F 1.8 mm, Allen Key, Ankylos	[]

	47.525.50 Driver Guide	[]
	47.525.55 Ratchet with demountable handle	[]
	47.803.02 Torque Wrench (optional)	[]

ALL THE TRADEMARKS ARE THE PROPERTY OF THE RESPECTIVE COMPANIES.
We assume no liability for deviations due to tolerances of implants.



HELMUT ZEPP
MEDIZINTECHNIK GMBH

Tel.: +49 (0) 7464/9888 - 0



HELMUT ZEPP
MEDIZINTECHNIK GMBH

www.zepf-dental.com