

ZEPF bonemill () () ()

Bone Mill




47.954.05 ZEPF Bone Mill Set

with helicoidal milling part

47.954.35

NEW!



NEW!



47.954.01 ZEPF Bone Mill

ZEPF bonemill () () ()

Bone Mill to crush autologous bone

To correct bone defects, often an interior bone removing on another place is necessary. Bigger bone pieces or boneblocks can be removed by using trephines.

The ZEPF Bone Mill allows to crush bone in order to produce the greatest possible volume of bone graft.

The grainy consistence of the produced bone graft guarantees an optimal adaption on the bone.

The extension bar which can be adapted on the rotary handle facilitates an optimized power transmission and torque.

The new Bone Mill with helicoidal milling part makes milling easier.

► **This new milling part can be ordered and subsequently be mounted in an older type of Bone Mill by the user himself.**

The ZEPF Bone Mill combines following advantages:

- easy handling
- quick assembly / disassembly without additional tools
- no loss of bone material (even in the case of small quantities)
- easy cleaning
- sterilization possible in washtray (autoclave)
- storage of sterile goods in container






HELMUT ZEPF
MEDIZINTECHNIK GMBH

NEWS

03/13


Distributed by:



HELMUT ZEPP
MEDIZINTECHNIK GMBH



Optional Instruments - Trephines

Article number	Ø inside	Ø outside	Teeth	Body length	Grading	Delivery quantity	Order quantity
08.910.01	1.7	2.3	7	22 mm	7/10/13/16	1 piece	[]
08.910.02	2.3	2.8	7	22 mm	7/10/13/16	1 piece	[]
08.910.03	2.8	3.3	9	22 mm	7/10/13/16	1 piece	[]
08.910.04	3.3	3.8	9	22 mm	7/10/13/16	1 piece	[]
08.910.05	4.0	4.5	11	22 mm	7/10/13/16	1 piece	[]
08.910.06	4.3	4.8	11	22 mm	7/10/13/16	1 piece	[]
08.910.07	4.8	5.8	9	22 mm	7/10/13/16	1 piece	[]
08.910.13	5.0	6.0	11	22 mm	7/10/13/16	1 piece	[]
08.910.08	6.0	7.0	11	22 mm	7/10/13/16	1 piece	[]
08.910.09	7.0	9.0	17	22 mm	7/10/13/16	1 piece	[]
08.910.10	9.0	9.0	17	22 mm	7/10/13/16	1 piece	[]
08.910.11	9.0	10.0	17	22 mm	7/10/13/16	1 piece	[]
08.910.12	10.0	11.0	19	22 mm	7/10/13/16	1 piece	[]
85.070.01	 Teflon® bur stand for 6 burs, 60° click into place, 85 x 50 mm					1 piece (without burs)	[]

Article description

Order quantity

47.954.05 ZEPP Bone Mill Set, consisting of:
Washtray incl. cleaning instrument
and Kirsch Sinus Plugger

[]

85.185.55 Washtray System
with Teflon® storage / holding system and lid
for ZEPP Bone Mill (without content)

[]

47.954.01 ZEPP Bone Mill,
bone mill to crush autologous bone,
as illustrated, incl. milling part

[]

47.954.35 Helicoidal Milling Part for Bone Mill 47.954.01
and to be mounted in 47.954.00 / 47.954.01

[]

47.954.55 Extension bar for the rotary handle 47.954.50,
to increase the lever action

[]

47.530.00 Cleaning instrument ZEPP-Line
to collect autologous bone

[]

41.868.07 Plugger acc. to Kirsch,
Titanium ZEPP-Line Handle, Sinus 7, double-ended,
17.5 cm, Spoon 8x10

[]

Order by fax to +49 (0) 7464 / 98 88 88 or online at: www.zepf-dental.com

Client Number

Company / Practice

Street, Number

Postal Code, City

Date

Source of Supply / Depot

Phone / Fax

Signature of the Client