



Rack for Endo Instruments



Endo Rack and Measuring Device for up to 18 Files or other Endo Instruments

Endo Rack, stainless steel, containing 2 cups and 18 drill holes.

The depth of the endo instruments can be adjusted exactly by means of the scale.

Blasted version, dimensions 87 x 82 x 20 mm.

NEWS

11/12



Article description

99.090.00 Endo Rack for up to 18 Endo Instruments

Order quantity

[]



HELMUT ZEPF
MEDIZINTECHNIK GMBH

78606 Seitingen-Oberflacht, Germany | Tel.: +49 (0) 7464 / 98 88 -0 | Fax: +49 (0) 7464 / 98 88 -88 | E-Mail: info@zepf-dental.com | Internet: www.zepf-dental.com

Order by fax to +49 (0) 7464 / 98 88 88 or online at: www.zepf-dental.com

Customer Number

Company / Practice

Street, Number

Postal Code, City

Date

Source of Supply / Depot

Phone / Fax

Signature of the Client